

Max. Marks : 19 Marks

Time : 19 Minutes

Four compounds, all colourless liquids, are

- butan-2-ol
- butanal
- butanone
- 2-methylpropan-2-ol

Two of these compounds can be identified using different test-tube reactions.


Describe these **two** test-tube reactions by giving reagents and observations in each case. Suggest how the results of a spectroscopic technique could be used to distinguish between the **other** two compounds.

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

(Total 6 marks)

Q2.


Cyclohexene (boiling point = 83 °C) can be prepared by the dehydration of cyclohexanol (boiling point = 161 °C) using concentrated phosphoric acid.



A student prepared cyclohexene by placing 10 cm³ of cyclohexanol (density = 0.96 g cm⁻³) into a round-bottomed flask.

3 cm³ of concentrated phosphoric acid were then carefully added to the flask.

The student added a few anti-bumping granules and set up the apparatus shown in the diagram.



- The student heated the mixture and collected the liquid that distilled at temperatures below 100 °C
- The distillate was poured into a separating funnel and washed by shaking with sodium carbonate solution.
- Periodically, the separating funnel was inverted and the tap opened.
- The aqueous layer was discarded and the final organic product was dried using anhydrous calcium chloride.
- After the product was dried, the drying agent was removed by filtration under reduced pressure.

(a) The student collected 5.97 g of cyclohexene in the experiment.

Calculate the percentage yield of cyclohexene.

(3)

- (b) Describe a test-tube reaction, on the product, to show that the cyclohexanol had been dehydrated.

State what you would observe. _____

(2)

- (c) Suggest why sodium carbonate solution was used to wash the distillate.

(1)

- (d) Explain why it is important to open the tap of the separating funnel periodically.

(1)

- (e) Give a property of anhydrous calcium chloride, other than its ability to absorb water, that makes it suitable as a drying agent in this preparation.

(1)

- (f) Describe the apparatus used to remove the drying agent by filtration under reduced pressure. Your description of the apparatus can be either a labelled diagram or a description in words.

(2)

- (g) A sample of cyclohexene has been contaminated with cyclohexanol. The cyclohexene can be separated from the cyclohexanol by column chromatography. Silica gel is used as the stationary phase and hexane as the mobile phase.

Explain why cyclohexene has a shorter retention time than cyclohexanol.

(2)

- (h) Explain how an infrared spectrum would confirm that the cyclohexene obtained from the chromatography column did not contain any cyclohexanol.

(1)

(Total 13 marks)